

# BCL2A1 Antibody



PACO07993

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

BCL2A1

**Uniprot**

Q16548

**Synonyms:**

BCL2-related protein A1; BCL2A1; ACC-1; ACC-2; BCL2L5; BFL1; GRS; HBPA1 ;

**Immunogen:**

Human BCL2A1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A