

BCL2L13 Antibody



PACO07995

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

BCL2L13

Uniprot

Q9BXK5

Synonyms:

BCL2-like 13 (apoptosis facilitator); BCL2L13; BCL-RAMBO; MIL1 ;

Immunogen:

Human BCL2L13.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A