

BCOR Antibody



PACO07999

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

BCOR

Uniprot

Q6W2J9

Synonyms:

BCL6 co-repressor; BCOR; ANOP2; FLJ20285; FLJ38041; KIAA1575; MAA2; MCOPS2; MGC131961; MGC71031 ;

Immunogen:

Human BCOR.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A