

BEST3 Antibody



PACO08005

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

BEST3

Uniprot

Q8N1M1

Synonyms:

bestrophin 3; BEST3; MGC13168; MGC40411; VMD2L3 ;

Immunogen:

Human BEST3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A