

BHLHE40 Antibody



PACO08009

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

BHLHE40

Uniprot

O14503

Synonyms:

basic helix-loop-helix family, member e40; BHLHE40; BHLHB2; DEC1; FLJ99214; SHARP-2; STRA13; Stra14 ; SHARP2

Immunogen:

Human BHLHE40.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A