

# BLVRB Antibody



PACO08024

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

BLVRB

**Uniprot**

P30043

**Synonyms:**

biliverdin reductase B (flavin reductase (NADPH)); BLVRB; BVRB; FLR; MGC117413; SDR43U1 ;

**Immunogen:**

Human BLVRB.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A