

BNIPL Antibody



PACO08033

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

BNIPL

Uniprot

Q7Z465

Synonyms:

BCL2/adenovirus E1B 19kD interacting protein like; BNIPL; BNIP-S; BNIPL-1; BNIPL-2; BNIPL1; BNIPL2; PP753 ;

Immunogen:

Human BNIPL

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A