

CELF5 Antibody



PACO08049

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CELF5

Uniprot

Q8N6W0

Synonyms:

bruno-like 5, RNA binding protein (Drosophila); BRUNOL5; BRUNOL-5; CELF5 ;

Immunogen:

Human BRUNOL5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A