

PACO08051

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

BSG

**Uniprot**

P35613

**Synonyms:**

basigin (Ok blood group); BSG;5F7; CD147; EMMPRIN; M6; OK; TCSF ;

**Immunogen:**

Human BSG.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A