

BTBD1 Antibody



PACO08053

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

BTBD1

Uniprot

Q9H0C5

Synonyms:

BTB (POZ) domain containing 1; BTBD1; C15orf1; NS5ATP8 ;

Immunogen:

Human BTBD1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A