

# BTN2A2 Antibody



PACO08061

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

BTN2A2

**Uniprot**

Q8WVV5

**Synonyms:**

butyrophilin, subfamily 2, member A2; BTN2A2; BT2.2; BTF2; FLJ41908 ;

**Immunogen:**

Human BTN2A2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A