

BUB1B Antibody



PACO08063

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

BUB1B

Uniprot

O60566

Synonyms:

budding uninhibited by benzimidazoles 1 homolog beta (yeast); BUB1B; BUB1beta; BUBR1; Bub1A; MAD3L; SSK1; hBUBR1 ;

Immunogen:

Human BUB1B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A