

# BUD31 Antibody



PACO08065

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

BUD31

**Uniprot**

P41223

**Synonyms:**

BUD31 homolog (S. cerevisiae); BUD31; EDG-2; EDG2; G10; MGC111202; YCR063W; fSAP17 ;

**Immunogen:**

Human BUD31.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A