

C11orf49 Antibody



PACO08071

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

C11orf49

Uniprot

Q9H6J7

Synonyms:

chromosome 11 open reading frame 49; C11orf49; FLJ22210; MGC4707 ;

Immunogen:

Human C11orf49.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A