

C11orf70 Antibody



PACO08074

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

C11orf70

Uniprot

Q9BRQ4

Synonyms:

chromosome 11 open reading frame 70; C11orf70; MGC13040 ;

Immunogen:

Human C11orf70.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A