

C11orf73 Antibody



PACO08075

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

HIKESHI

Uniprot

Q53FT3

Synonyms:

chromosome 11 open reading frame 73; C11orf73; FLJ43020; HSPC138; HSPC179; L7RN6;

Immunogen:

Human C11orf73.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A