

KIAA0226L Antibody



PACO08077

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

RUBCNL

Uniprot

Q9H714

Synonyms:

chromosome 13 open reading frame 18; C13orf18; FLJ21562; FLJ43762 ;

Immunogen:

Human C13orf18.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A