

VIPAS39 Antibody



PACO08078

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

VIPAS39

Uniprot

Q9H9C1

Synonyms:

chromosome 14 open reading frame 133; C14orf133; FLJ12707 ; VIPAR

Immunogen:

Human C14orf133.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A