

HAUS4 Antibody



PACO08080

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

HAUS4

Uniprot

Q9H6D7

Synonyms:

chromosome 14 open reading frame 94; C14orf94; FLJ20424 ;

Immunogen:

Human C14orf94.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A