

FAM195A Antibody



PACO08081

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

MCRIP2

Uniprot

Q9BUT9

Synonyms:

chromosome 16 open reading frame 14; C16orf14; MGC104301; MGC15416; c349E10.1
;; FAM195A

Immunogen:

Human C16orf14.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A