

HID1 Antibody



PACO08083

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

HID1

Uniprot

Q8IV36

Synonyms:

chromosome 17 open reading frame 28; C17orf28; DMC1; FLJ43526 ;

Immunogen:

Human C17orf28.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A