

PACO08091

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SARG

Uniprot

Q9BW04

Synonyms:

chromosome 1 open reading frame 116; C1orf116; DKFZp666H2010; FLJ36507; MGC2742; MGC4309; SARG ;

Immunogen:

Human C1orf116.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A