

PPDPF Antibody



PACO08101

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

PPDPF

Uniprot

Q9H3Y8

Synonyms:

chromosome 20 open reading frame 149; C20orf149; FLJ21046; MGC2479; dJ697K14.9 ;
PPDPF, EXPDF

Immunogen:

Human C20orf149.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A