

# C9orf156 Antibody



PACO08117

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TRMO

**Uniprot**

Q9BU70

**Synonyms:**

chromosome 9 open reading frame 156; C9orf156; HSPC219; NAP1 ;

**Immunogen:**

Human C9orf156.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A