

CA11 Antibody



PACO08121

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CA11

Uniprot

O75493

Synonyms:

carbonic anhydrase XI; CA11; CARP2; CARPX1 ;

Immunogen:

Human CA11.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A