

# CAB39L Antibody



PACO08131

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:****Protein Background:****Gene ID:**

CAB39L

**Uniprot**

Q9H9S4

**Synonyms:**

calcium binding protein 39-like; CAB39L; FLJ12577; MO25-BETA; MO2L; RP11-103J18.3; bA103J18.3 ;

**Immunogen:**

Human CAB39L.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A