

PACO08133

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CABYR

Uniprot

O75952

Synonyms:

calcium binding tyrosine-(Y)-phosphorylation regulated; CABYR; CBP86; FSP-2; FSP2; MGC9117 ;

Immunogen:

Human CABYR.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A