

CACNG3 Antibody



PACO08136

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CACNG3

Uniprot

O60359

Synonyms:

calcium channel, voltage-dependent, gamma subunit 3; CACNG3; Cacng2 ;

Immunogen:

Human CACNG3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A