

CAP1 Antibody



PACO08165

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CAP1

Uniprot

Q01518

Synonyms:

CAP, adenylate cyclase-associated protein 1 (yeast); CAP1; CAP; CAP1-PEN ;

Immunogen:

Human CAP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A