

CAPN9 Antibody



PACO08175

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

CAPN9

Uniprot

O14815

Synonyms:

calpain 9; CAPN9; GC36; nCL-4 ;

Immunogen:

Human CAPN9.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A