

CAPS2 Antibody



PACO08180

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CAPS2

Uniprot

Q9BXY5

Synonyms:

calcyphosine 2; CAPS2; FLJ34520; UG0636c06 ;

Immunogen:

Human CAPS2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A