

CASC3 Antibody



PACO08191

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:**Protein Background:****Gene ID:**

CASC3

Uniprot

O15234

Synonyms:

cancer susceptibility candidate 3; CASC3; BTZ; MLN51 ;

Immunogen:

Human CASC3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A