

PACO08192

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CASK

Uniprot

O14936

Synonyms:

calcium/calmodulin-dependent serine protein kinase (MAGUK family); CASK; CAGH39; CMG; FGS4; FLJ22219; FLJ31914; LIN2; MICPCH; TNRC8 ;

Immunogen:

Human CASK.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A