

# CC2D1B Antibody



PACO08219

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CC2D1B

**Uniprot**

Q5T0F9

**Synonyms:**

coiled-coil and C2 domain containing 1B; CC2D1B; KIAA1836; RP11-155O18.2 ;

**Immunogen:**

Human CC2D1B.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A