

# CCDC115 Antibody



PACO08221

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CCDC115

**Uniprot**

Q96NT0

**Synonyms:**

coiled-coil domain containing 115; CCDC115; FLJ30131; MGC12981; ccp1 ;

**Immunogen:**

Human CCDC115.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A