

CCDC158 Antibody



PACO08223

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

CCDC158

Uniprot

Q5M9N0

Synonyms:

coiled-coil domain containing 158; CCDC158; FLJ25770; MGC35086 ;

Immunogen:

Human CCDC158.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A