

# CCDC28A Antibody



PACO08226

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CCDC28A

**Uniprot**

Q8IWP9

**Synonyms:**

coiled-coil domain containing 28A; CCDC28A; C6orf80; CCRL1AP; DKFZp586D0623; MGC131913 ;

**Immunogen:**

Human CCDC28A.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A