

CCDC28B Antibody



PACO08227

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CCDC28B

Uniprot

Q9BUN5

Synonyms:

coiled-coil domain containing 28B; CCDC28B; MGC1203; MGC16441; RP4-622L5.5 ;

Immunogen:

Human CCDC28B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A