

COA3 Antibody



PACO08232

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

COA3

Uniprot

Q9Y2R0

Synonyms:

coiled-coil domain containing 56; CCDC56; FLJ50764; HSPC009 ;

Immunogen:

Human CCDC56.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A