

CCL11 Antibody



PACO08239

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CCL11

Uniprot

P51671

Synonyms:

chemokine (C-C motif) ligand 11; CCL11; MGC22554; SCYA11 ;

Immunogen:

Human CCL11.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A