

# CCL28 Antibody



PACO08241

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CCL28

**Uniprot**

Q9NRJ3

**Synonyms:**

chemokine (C-C motif) ligand 28; CCL28; CCK1; MEC; MGC71902; SCYA28 ;

**Immunogen:**

Human CCL28.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A