

CCNL1 Antibody



PACO08255

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CCNL1

Uniprot

Q9UK58

Synonyms:

cyclin L1; CCNL1; BM-001; PRO1073; ania-6a ;

Immunogen:

Human CCNL1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A