

CD1A Antibody



PACO08271

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:**Protein Background:****Gene ID:**

CD1A

Uniprot

P06126

Synonyms:

CD1a molecule; CD1A; CD1; FCB6; HTA1; R4; T6 ;

Immunogen:

Human CD1A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A