

CD248 Antibody



PACO08277

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

CD248

Uniprot

Q9HCU0

Synonyms:

CD248 molecule, endosialin; CD248; CD164L1; MGC119478; MGC119479; TEM1 ;

Immunogen:

Human CD248.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A