

CD3E Antibody



PACO08287

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CD3E

Uniprot

P07766

Synonyms:

CD3e molecule, epsilon (CD3-TCR complex); CD3E; FLJ18683; T3E; TCRE ;

Immunogen:

Human CD3E.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A