

CD82 Antibody



PACO08301

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CD82

Uniprot

P27701

Synonyms:

CD82 molecule; CD82;4F9; C33; GR15; IA4; KAI1; R2; SAR2; ST6; TSPAN27 ;

Immunogen:

Human CD82.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A