

# CDC25B Antibody



PACO08311

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

CDC25B

**Uniprot**

P30305

**Synonyms:**

cell division cycle 25 homolog B (S. pombe);

**Immunogen:**

Human CDC25B.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A