

# CDC37L1 Antibody



PACO08317

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CDC37L1

**Uniprot**

Q7L3B6

**Synonyms:**

cell division cycle 37 homolog (S. cerevisiae)-like 1; CDC37L1; CDC37B; FLJ20639; HARC ;

**Immunogen:**

Human CDC37L1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A