

CDC42EP2 Antibody



PACO08319

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CDC42EP2

Uniprot

O14613

Synonyms:

CDC42 effector protein (Rho GTPase binding) 2; CDC42EP2; BORG1; CEP2 ;

Immunogen:

Human CDC42EP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A