

CDC73 Antibody



PACO08326

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

CDC73

Uniprot

Q6P1J9

Synonyms:

cell division cycle 73, Paf1/RNA polymerase II complex component, homolog (S. cerevisiae); CDC73; C1orf28; FLJ23316; HPT-JT; HRPT2 ;

Immunogen:

Human CDC73.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A